



**Memo No.** 2018-046  
**Date:** 22-Oct-2018  
**Memo To:** Clients  
**Re:** Test Revisions

---

***Antithrombin III Antigen, Plasma (Citrate)*** Revised method, reference range and units. Refer to website for details. This method is a quantitative assay with enhanced sensitivity, specificity and reproducibility.

---

***Ammonia, Plasma (EDTA)***

The reference ranges listed on the website have been revised:

Adult: 0 - 48 µmol/L (previously listed as "Less than 48")

Neonate: 0 – 73 µmol/L (previously listed as "Less than 73")

---

***C-Peptide Stimulation***

Assay discontinued.

---

***Citrate, Urine and Oxalate, Urine***

Specimen handling instructions have been revised to request that separate specimens be submitted for each of these tests.

---

***α-Fetoprotein, Serum***

Pediatric reference ranges have been added. Biotin supplements can affect results and an advisory has been added. Refer to the website for details.

---

***Mucopolysaccharides, Urine & Oligosaccharides Screen, Urine***

Collection instructions have been revised for both tests to remove instruction to avoid first morning collection. A random collection at any time is acceptable.

---

If you have further questions, please contact Client Care at (416) 422-3000 Ext. 300 or [info@ICLabs.ca](mailto:info@ICLabs.ca)

Dr. Shashank Tilak, Ph.D  
Laboratory Director  
(416) 422-3000 Ext. 221  
[ShashTilak@ICLabs.ca](mailto:ShashTilak@ICLabs.ca)

**Want to receive updates by e-mail? Please contact Client Care [info@ICLabs.ca](mailto:info@ICLabs.ca)**