



Memo No. 2018-026
Date: 29-Jun-2018
Memo To: Clients
Re: Test Revision – Dihydrotestosterone

Please note that reference values will be revised effective 03-Jul-2018 with the addition of pediatric values and a small revision for the lower end of the range for post-menopausal women.

There is no change in method or test site.

The conversion factor listed on the ICL website has been corrected to pmol/L = ng/dL x 34.4

Revised Dihydrotestosterone Reference Values (Effective 3-Jul-2018)

Male:

< 8 years: ≤ 632 pmol/L

8 – 18 years: ≤ 3367 pmol/L

adult: 860 – 3406 pmol/L

Female:

< 8 years: ≤ 632 pmol/L

8 – 18 years: ≤ 1916 pmol/L

premenopausal adult: 83 – 1266 pmol/L

postmenopausal: 34 – 623 pmol/L

If you have further questions, please contact Client Care at (416) 422-3000 Ext. 300 or info@ICLabs.ca or you may contact me directly.

Shashank Tilak PhD, DCC
Laboratory & Scientific Director
(416) 422-3000 Ext. 221
ShashTilak@ICLabs.ca

Want to receive updates by e-mail? Please contact Client Care info@ICLabs.ca

ICL is the only Not-For-Profit dedicated broker of lab diagnostic tests in Canada

57 Gervais Drive, North York, ON M3C 1Z2

(416) 422-3000 Ext. 300 info@ICLabs.ca